




CYP1B1 Antibody

Product Code	CSB-PA273340
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16678
Immunogen	Synthesized peptide derived from C-terminal of Human CP1B1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CP1B1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	aryl hydrocarbon hydroxylase; CP1B; CP1B1; CYP1B1; CYP1B1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP1B1

Image



Western blot analysis of extracts from Jurkat cells, using CP1B1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.