




PTGIS Antibody

Product Code	CSB-PA955037
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16647
Immunogen	Synthesized peptide derived from Internal of Human PTGIS.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total PTGIS protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CYP8; CYP8A1; cytochrome P450; family 8; subfamily A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTGIS

Image



Western blot analysis of extracts from K562 cells, using PTGIS antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.