




TBXAS1 Antibody

Product Code	CSB-PA589726
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24557
Immunogen	Synthesized peptide derived from C-terminal of Human THAS.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total THAS protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CYP5; CYP5A1; Cytochrome P450 5A1; cytochrome P450; family 5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TBXAS1

Image



Western blot analysis of extracts from Jurkat cells, using THAS antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.