




CYP3A7 Antibody

Product Code	CSB-PA929315
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4D288
Immunogen	Synthesized peptide derived from internal of Human Cytochrome P450 3A7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 3A7 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cytochrome P450; family 3; subfamily A; polypeptide 7; CYP3A7
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP3A7

Image



Western blot analysis of extracts from LOVO cells, using Cytochrome P450 3A7 antibody.