




CYP39A1 Antibody

Product Code	CSB-PA932307
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYL5
Immunogen	Synthesized peptide derived from Internal of Human CYP39A1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CYP39A1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	24-hydroxycholesterol 7-alpha-hydroxylase; CP39A; CYP39A1; Cytochrome P450 39A1; cytochrome P450
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP39A1

Image



Western blot analysis of extracts from 293/HeLa cells, using CYP39A1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.