





CYP8B1 Antibody

Product Code	CSB-PA587317
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNU6
Immunogen	Synthesized peptide derived from internal of Human Cytochrome P450 8B1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 8B1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cytochrome P450; family 8; subfamily B; polypeptide 1; CYP8B1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP8B1

Image



Western blot analysis of extracts from 293 cells, using Cytochrome P450 8B1 antibody.



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue, using Cytochrome P450 8B1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.