





BAG3 Antibody

Product Code	CSB-PA644586
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95817
Immunogen	Synthesized peptide derived from internal of Human BAG3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total BAG3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BAG-3; BAG-family molecular chaperone regulator-3; BCL-2 binding athanogene- 3; Bcl-2-binding protein Bis; BIS
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BAG3

Image



Western blot analysis of extracts from 293 cells, using BAG3 antibody.



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue, using BAG3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.