




BAG4 Antibody

Product Code	CSB-PA555254
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95429
Immunogen	Synthesized peptide derived from C-terminal of Human BAG4.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous levels of total BAG4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BAG family molecular chaperone regulator 4; BAG-4; BAG-family molecular chaperone regulator-4; BAG4; Bcl-2-associated athanogene 4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BAG4

Image



Western blot analysis of extracts from K562 cells, using BAG4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.