




BAG4 Antibody

Product Code	CSB-PA555254
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95429
Immunogen	Synthesized peptide derived from C-terminal of Human BAG4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total BAG4protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BAG family molecular chaperone regulator 4; BAG-4; BAG-family molecular chaperone regulator-4; BAG4; Bcl-2-associated athanogene 4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BAG4

Image



Western blot analysis of extracts from K562 cells, using BAG4 antibody.