





STK17B Antibody

Product Code	CSB-PA574347
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94768
Immunogen	Synthesized peptide derived from internal of Human STK17B.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total STK17B protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DAP kinase-related apoptosis-inducing protein kinase 2; DAP kinase-related apoptotic kinase 2; Death-associated protein kinase-related apoptosis inducing protein kinase; kinase DRAK2; Serine/threonine kinase 17B
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STK17B

Image



Western blot analysis of extracts from mouse brain cells, using STK17B antibody.



Immunohistochemistry analysis of paraffin-embedded human tonsil tissue, using STK17B antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.