




EPHB1 Antibody

Product Code	CSB-PA073597
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54762
Immunogen	Synthesized peptide derived from C-terminal of Human EPHB1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total EPHB1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ELK; EPB1; Eph receptor B1; Ephrin type-B receptor 1 precursor; EPHT2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EPHB1

Image



Western blot analysis of extracts from COS7 cells, using EPHB1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.