




STK10 Antibody

Product Code	CSB-PA889950
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94804
Immunogen	Synthesized peptide derived from internal of Human STK10.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total STK10 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	kinase GEK2; LOK; Lymphocyte-oriented kinase; Serine/threonine protein kinase 10; STK10
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STK10

Image



Western blot analysis of extracts from HeLa cells and HuvEc cells, using STK10 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.