




MERTK Antibody

Product Code	CSB-PA253602
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12866
Immunogen	Synthesized peptide derived from N-terminal of Human MERTK.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MERTK protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	C-mer; kinase Mer; MERK; MERTK; proto-oncogene tyrosine-protein kinase MER precursor
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MERTK

Image



Western blot analysis of extracts from K562 cells and Jurkat cells, using MERTK antibody.