




CDC42BPA Antibody

Product Code	CSB-PA019608
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VT25
Immunogen	Synthesized peptide derived from Internal of Human MRCKa.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:200,IF:1:100-1:200
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDC42BPA

Image



Western blot analysis of extracts from Jurkat/COLO205 cells, using MRCKa antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.