





PAK7 Antibody

Product Code	CSB-PA265580
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P286
Immunogen	Synthesized peptide derived from C-terminal of Human PAK7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PAK7 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	KIAA1264; p21-activated kinase 7; PAK 7; PAK-5; PAK-7
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PAK5

Image



Western blot analysis of extracts from rat brain cells, using PAK7 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using PAK7 antibody.