




PHKG1 Antibody

Product Code	CSB-PA046995
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16816
Immunogen	Synthesized peptide derived from internal of Human PHKG1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PHKG1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	KPBG; PHKG; Phosphorylase B kinase gamma catalytic chain skeletal muscle isoform; Phosphorylase kinase gamma subunit 1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHKG1

Image



Western blot analysis of extracts from HuvEc cells and HepG2 cells, using PHKG1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.