




CD79A Antibody

Product Code	CSB-PA165334
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11912
Immunogen	Synthesized peptide derived from internal of Human CD79A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CD79A protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	B-cell antigen receptor complex-associated protein alpha chain; Ig-alpha; MB-1 membrane glycoprotein; Surface IgM-associated protein; Membrane-bound immunoglobulin-associated protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD79A

Image



Western blot analysis of extracts from HuvEc cells, using CD79A antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.