




THY1 Antibody

Product Code	CSB-PA170582
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04216
Immunogen	Synthesized peptide derived from internal of Human THY1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total THY1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Thy-1 membrane glycoprotein; Thy-1 antigen; CDw90; CD90; THY1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	THY1

Image



Western blot analysis of extracts from HuvEc cells, using THY1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.