




LAMP1 Antibody

Product Code	CSB-PA616273
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11279
Immunogen	Synthesized peptide derived from internal of Human LAMP1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total LAMP1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Lysosome-associated membrane glycoprotein 1; Lysosome-associated membrane protein 1; LAMP-1; CD107 antigen-like family member A; CD107a
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LAMP1

Image



Western blot analysis of extracts from K562, Jurkat, HepG2, HuvEc, HeLa and COLO205 cells, using LAMP1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.