




TNFRSF10D Antibody

Product Code	CSB-PA031451
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBN6
Immunogen	Synthesized peptide derived from internal of Human TNFRSF10D.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TNFRSF10D protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 10D; Decoy receptor 2; DcR2; TNF-related apoptosis-inducing ligand receptor 4; TRAIL receptor 4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFRSF10D

Image



Western blot analysis of extracts from Jurkat cells, using TNFRSF10D antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.