




BTLA Antibody

Product Code	CSB-PA921165
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z6A9
Immunogen	Synthesized peptide derived from C-terminal of Human BTLA.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total BTLA protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	B and T lymphocyte associated; B and T lymphocyte attenuator; B- and T-lymphocyte attenuator; B- and T-lymphocyte-associated protein; BTLA
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BTLA

Image



Western blot analysis of extracts from K562/Jurkat cells, using BTLA antibody.