




FAM3D Antibody

Product Code	CSB-PA586397
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BQ1
Immunogen	Synthesized peptide derived from C-terminal of Human FAM3D.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total FAM3Dprotein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	cytokine-like protein EF-7; EF7; FAM3D; family with sequence similarity 3; member D
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FAM3D

Image



Western blot analysis of extracts from K562 cells, using FAM3D antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.