




TOLLIP Antibody

Product Code	CSB-PA571897
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0E2
Immunogen	Synthesized peptide derived from Internal of Human TOLLIP.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TOLLIPprotein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	adapter protein; FLJ33531; IL-1RAcPI; TOLIP; toll interacting protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TOLLIP

Image



Western blot analysis of extracts from HuvEc/Jurkat/293 cells, using TOLLIP antibody.