





AKT1 Antibody

Product Code	CSB-PA204523
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31749
Immunogen	Fusion protein corresponding to residues near the N terminal of Human V-akt murine thymoma viral oncogene homolog 1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total AKT1 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AKT1

Image



Immunohistochemical analysis of paraffin-embedded Human colon cancer tissue using at dilution 1/15.



Immunohistochemical analysis of paraffin-embedded Human tonsil tissue using at dilution 1/15.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.