




CYP11A1 Antibody

Product Code	CSB-PA914189
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05108
Immunogen	Synthesized peptide derived from C-terminal of Human Cytochrome P450 11A1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 11A1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cholesterol side-chain cleavage enzyme mitochondrial; CYPXIA1; Cholesterol desmolase; Cytochrome P450 11A1; Cytochrome P450(scc)
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP11A1

Image



Western blot analysis of extracts from HeLa cells and A549 cells, using Cytochrome P450 11A1 antibody.