





NEUROG3 Antibody

Product Code	CSB-PA956774
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4Z2
Immunogen	Full length fusion protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total NEUROG3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	ngn3; Atoh5; NGN-3; Math4B; bHLHa7
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NEUROG3

Image



Gel: 12%SDS-PAGE Lysates (from left to right):
 Mouse liver tissue Amount of lysate: 50µg per
 lane Primary antibody: 1/250 dilution Secondary
 antibody dilution: 1/8000 Exposure time: 40
 seconds



Immunohistochemical analysis of paraffin-
 embedded Human thyroid cancer tissue using at
 dilution 1/30.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.