




CTLA4 Antibody

Product Code	CSB-PA224600
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16410
Immunogen	Synthetic peptide corresponding to residues near the C terminal of Human Cytotoxic T-lymphocyte protein 4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CTLA4 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	CD; GSE; GRD4; ICOS; CD152; CTLA-4; IDDM12; CELIAC3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CTLA4

Image



Immunohistochemical analysis of paraffin-embedded Human cervical cancer tissue using at dilution 1/25.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.