




GSK3B Antibody

Product Code	CSB-PA939628
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49841
Immunogen	Synthetic peptide corresponding to a region derived from internal residues of Human glycogen synthase kinase 3 beta
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total GSK3B protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	GSK-3 beta; GSK3B; GSK3beta isoform
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GSK3B

Image



Gel: 10%SDS-PAGE Lysates (from left to right): Human small intestine stromal cancer tissue
 Amount of lysate: 30µg per lane Primary antibody: 1/100 dilution Secondary antibody dilution: 1/8000 Exposure time: 5 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.