




MC3R Antibody

Product Code	CSB-PA261897
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41968
Immunogen	Synthetic peptide corresponding to residues near the N terminal of Human melanocortin 3 receptor
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total MC3R protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	MC3; OB20; OQTL; BMIQ9; MC3-R
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MC3R

Image



Gel: 10%SDS-PAGE Lysates (from left to right):
 Mouse large intestine tissue Amount of lysate:
 30µg per lane Primary antibody: 1/800 dilution
 Secondary antibody dilution: 1/8000 Exposure
 time: 20 seconds