




# MC3R Antibody

<b>Product Code</b>	CSB-PA261897
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P41968
<b>Immunogen</b>	Synthetic peptide corresponding to residues near the N terminal of Human melanocortin 3 receptor
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous levels of total MC3R protein.
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:3000
<b>Form</b>	Rabbit IgG in pH7.3 PBS, 0.05% NaN <sub>3</sub> , 50% Glycerol.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Clonality</b>	Polyclonal
<b>Alias</b>	MC3; OB20; OQTL; BMIQ9; MC3-R
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MC3R

## Image



Gel: 10%SDS-PAGE Lysates (from left to right):  
 Mouse large intestine tissue Amount of lysate:  
 30µg per lane Primary antibody: 1/800 dilution  
 Secondary antibody dilution: 1/8000 Exposure  
 time: 20 seconds

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.