





# AKAP1 Antibody

|                            |                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA793193                                                                                                         |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                        |
| <b>Uniprot No.</b>         | Q92667                                                                                                               |
| <b>Immunogen</b>           | Synthetic peptide corresponding to a region derived from internal residues of Human A kinase (PRKA) anchor protein 1 |
| <b>Raised In</b>           | Rabbit                                                                                                               |
| <b>Species Reactivity</b>  | Human,Mouse,Rat                                                                                                      |
| <b>Specificity</b>         | The antibody detects endogenous levels of total AKAP1 protein.                                                       |
| <b>Tested Applications</b> | ELISA,IHC;IHC:1:50-1:100                                                                                             |
| <b>Form</b>                | Rabbit IgG in pH7.3 PBS, 0.05% NaN <sub>3</sub> , 50% Glycerol.                                                      |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                                                           |
| <b>Alias</b>               | AKAP; PRKA1; AKAP84; TDRD17; AKAP121; AKAP149; D-AKAP1; PPP1R43; SAKAP84                                             |
| <b>Product Type</b>        | Polyclonal Antibody                                                                                                  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                 |
| <b>Target Names</b>        | AKAP1                                                                                                                |

## Image



Immunohistochemical analysis of paraffin-embedded Human colon cancer tissue using at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using at dilution 1/40.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.