





SCN11A Antibody

Product Code	CSB-PA961122
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UI33
Immunogen	Synthetic peptide corresponding to a region derived from internal residues of Human sodium channel, voltage-gated, type XI, alpha subunit
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total SCN1A/2A/3A/4A/5A/8A/9A/10A/11A/12A protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	NaN; SNS-2; NAV1.9; SCN12A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SCN11A

Image



Immunohistochemical analysis of paraffin-embedded Human cervical cancer tissue using at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using at dilution 1/40.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.