





CACNA2D1 Antibody

Product Code	CSB-PA059719
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54289
Immunogen	Synthetic peptide corresponding to a region derived from internal residues of Human calcium channel, voltage-dependent, L type, alpha 1D subunit
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CACNA1 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	CACH3; CACN4; PASNA; SANDD; Cav1.3; CCHL1A2; CACNL1A2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CACNA2D1

Image



Immunohistochemical analysis of paraffin-embedded Human gastric cancer tissue using at dilution 1/30.



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using at dilution 1/30.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.