





KRT6C Antibody

Product Code	CSB-PA216033
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48668
Immunogen	Synthetic peptide corresponding to a region derived from internal residues of Human keratin 6C
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous levels of total KRT6C protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	K6E; KRT6E
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KRT6C

Image



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using.



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.