




MECP2 Antibody

Product Code	CSB-PA912786
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51608
Raised In	Rabbit
Species Reactivity	Human
Specificity	MeCP2 rabbit polyclonal antibody detects endogenous MeCP2 proteins.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	AUTSX3; DKFZp686A24160; MeCp-2 protein; MECP2; methyl CpG binding protein 2 (Rett syndrome)
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MECP2

Image



Western blot analysis of HeLa cell lysate, using diluted at 1:2,000 (lane 1) and 1:5,000 (lane 2).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.