




FH Antibody

Product Code	CSB-PA553013
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07954
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	Fumarase rabbit polyclonal antibody detects endogenous fumarase proteins.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	FH; Fumarase; fumarate hydratase; Fumarate hydratase, mitochondrial; FUMH; HLRCC; LRCC; MCL; MCUL1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FH

Image



Western blot analysis of Hela (lane 1), mouse brain (lane 2) and HepG2 (lane 3), using diluted at 1:2,000.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.