



AmCyan Mouse Polyclonal Antibody

Product Code	CSB-PA218792
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	The AmCyan Mouse Polyclonal Antibody detects AmCyan & AmCyan Tag protein.
Raised In	Mouse
Species Reactivity	Human Rat Mouse Monkey Dog Chicken Hamster Rabbit Pig Sheep
Tested Applications	WB, IHC, IP, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.