



# EBFP Mouse Polyclonal Antibody

<b>Product Code</b>	CSB-PA588045
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9ULX5
<b>Immunogen</b>	The EBFP Mouse Polyclonal Antibody detects EBFP & EBFP Tag protein.
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Tested Applications</b>	WB; WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Purification Method</b>	Antigen Affinity Purified
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Target Names</b>	RNF112