



EBFP Mouse Polyclonal Antibody

Product Code	CSB-PA588045
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULX5
Immunogen	The EBFP Mouse Polyclonal Antibody detects EBFP & EBFP Tag protein.
Raised In	Mouse
Species Reactivity	N/A
Tested Applications	WB; WB:1:500-1:2000
Form	Liquid
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Target Names	RNF112
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.