




# ECFP Mouse Polyclonal Antibody

<b>Product Code</b>	CSB-PA193774
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Raised In</b>	Mouse
<b>Specificity</b>	The ECFP Antibody detects ECFP & CFP tag protein.
<b>Tested Applications</b>	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
<b>Form</b>	Liquid
<b>Purification Method</b>	Antigen Affinity Purified
<b>Clonality</b>	Polyclonal

## Image



Western blot analysis of ECFP Recombinant protein using diluted at 1) 1:5000 2) 1:10000

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.