




ECFP Mouse Polyclonal Antibody

Product Code	CSB-PA193774
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Mouse
Specificity	The ECFP Antibody detects ECFP & CFP tag protein.
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal

Image



Western blot analysis of ECFP Recombinant protein using diluted at 1) 1:5000 2) 1:10000