



RFP Mouse Polyclonal Antibody

Product Code	CSB-PA978335
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Mouse
Species Reactivity	N/A
Specificity	The RFP tag antibody can recognize RFP tag fusion proteins.
Tested Applications	WB, IHC, IP, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, ELISA:1:20000
Form	Liquid
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Target Names	TRIM27