




GFP Mouse Polyclonal Antibody

Product Code	CSB-PA155763
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Mouse
Tested Applications	WB, IHC, IP, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Product Type	Polyclonal Antibody

Image



Western blot analysis of GFP Recombinant protein using diluted at 1) 1:5000 2) 1:10000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.