




# Phospho-GAB2 (Tyr452) Antibody

<b>Product Code</b>	CSB-PA781004
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UQC2
<b>Immunogen</b>	Peptide sequence around phosphorylation site of tyrosine 452 (D-N-Y(p)-V-P) derived from Human Gab2.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous levels of Gab2 only when phosphorylated at Tyr452.
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:3000
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
<b>Clonality</b>	Polyclonal
<b>Alias</b>	GRB2-associated-binding protein 2 GRB2-associated binder 2 Growth factor receptor bound protein 2-associated protein 2 pp100 GAB2 KIAA0571
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	GAB2

## Image



Western blot analysis of extracts from JK, using Gab2 (Phospho-Tyr452) antibody. The lane on the right is treated with the synthesized peptide.

<b>Product Modify</b>	Phospho-Tyr452
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.