





Phospho-UBTF (Ser484) Antibody

Product Code	CSB-PA072052
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17480
Immunogen	Peptide sequence around phosphorylation site of serine 484 (P-E-S(p)-P-K) derived from Human UBF.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of UBF only when phosphorylated at serine 484.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	Autoantigen NOR-90; Nucleolar transcription factor 1; TCFUBF; UBF-1; UBF1; UBTF; Upstream binding factor 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBTF

Image



Western blot analysis of extracts from HeLa cells, treated with calyculinA (50ng/ml, 30mins), using UBF (Phospho-Ser484) antibody. The lane on the right is treated with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using UBF (Phospho-Ser484) antibody. The picture on the right is treated with the synthesized peptide.

Product Modify

Phospho-Ser484

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.