




CD7 Antibody

Product Code	CSB-PA004953ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09564
Immunogen	Recombinant Human T-cell antigen CD7 protein (26-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell antigen CD7 (GP40) (T-cell leukemia antigen) (T-cell surface antigen Leu-9) (TP41) (CD antigen CD7), CD7
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD7

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA004953ESR1HU at dilution of 1:100