





AKR7A3 Antibody

Product Code	CSB-PA001551ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95154
Immunogen	Recombinant Human Aflatoxin B1 aldehyde reductase member 3 protein (1-331AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aflatoxin B1 aldehyde reductase member 3 (EC 1.-.-.) (AFB1 aldehyde reductase 2) (AFB1-AR 2), AKR7A3, AFAR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AKR7A3

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA001551ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA001551ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.