





CCK Antibody

Product Code	CSB-PA004771ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06307
Immunogen	Recombinant Human Cholecystokinin protein (21-115AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cholecystokinin (CCK) [Cleaved into: Cholecystokinin-58 (CCK58); Cholecystokinin-58 desnonopeptide ((1-49)-CCK58); Cholecystokinin-39 (CCK39); Cholecystokinin-33 (CCK33); Cholecystokinin-25 (CCK25); Cholecystokinin-18 (CCK18); Cholecystokinin-12 (CCK12); Cholecystokinin-8 (CCK8); Cholecystokinin-7 (CCK7); Cholecystokinin-5 (CCK5)], CCK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CCK

Image



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA004771ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA004771ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.