




CHMP1B Antibody

Product Code	CSB-PA748726ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7LBR1
Immunogen	Recombinant Human Charged multivesicular body protein 1b protein (1-199AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Charged multivesicular body protein 1b (CHMP1.5) (Chromatin-modifying protein 1b) (CHMP1b) (Vacuolar protein sorting-associated protein 46-2) (Vps46-2) (hVps46-2), CHMP1B, C18orf2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CHMP1B

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA748726ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.