





AGGF1 Antibody

Product Code	CSB-PA001444ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N302
Immunogen	Recombinant Human Angiogenic factor with G patch and FHA domains 1 protein (1-70AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiogenic factor with G patch and FHA domains 1 (Angiogenic factor VG5Q) (hVG5Q) (G patch domain-containing protein 7) (Vasculogenesis gene on 5q protein), AGGF1, GPATC7 GPATCH7 VG5Q
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	AGGF1

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA001444ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA001444ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.